

# AD302L

Air Drying

Part number: 40-1220




## Description

The AD302L is a useful attachment for any laboratory's external source of air, removing moisture and ensuring the delivery of dry air. Moisture content in house air can contaminate lab instrumentation, impacting on results and in some cases damaging the instrument itself. Compact in size, wall-mountable and light, the AD302L requires little space and its few moving parts mean minimal maintenance is required.



## Generator Details

- **Gas types:** Air
- **Integrated compressor:** No



## Key Features

- Produces clean, dry air, ensures optimum supply for gas generators and instruments
- Dryer performs 24/7 if needed
- Dryer is wall-mountable for small labs
- Few moving parts means little maintenance required and ensures long life of the generator
- Completely silent in operation
- Minimal set-up required
- 12 month comprehensive on-site warranty

## Technical Specifications

### AD302L - Air Dryer

|                      |                          |
|----------------------|--------------------------|
| Gas Type             | Air                      |
| Max Gas Flow:        | 302L/min                 |
| Gas Outlets Fitting: | 1x 1/4" BSP Female       |
| Min Inlet Pressure:  | 100psi/6.9bar            |
| Max Inlet Pressure:  | 120psi/8.3bar            |
| Start Up Time:       | 90 mins                  |
| Power Consumption:   | 8 - 19.2 W               |
| Voltage:             | 100 - 240 V              |
| Frequency:           | 50 / 60 Hz               |
| Current:             | 0.8 A                    |
| Accreditations:      | CE                       |
| Size (HxWxD) mm      | 910 x 160 x 471 mm       |
| Size (HxWxD) Inches  | 35.8 x 6.3 x 18.5 inches |
| Generator Weight     | 20.5kg / 45.1lbs         |

# [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

