

# TOC 1500

Zero Air Co2 Free

Part number: 42-1030



## Description

Designed to provide a reliable air supply for TOC analyzers, the TOC1500HP is a robust system used in labs to produce clean, dry CO<sub>2</sub> and CH<sub>4</sub> free




## Application Types

TOC Analyzers

## Generator Details

- **Gas types:** Air
- **Integrated compressor:** No



## Key Features

- Produces clean, dry, CO<sub>2</sub> and CH<sub>4</sub>-free air
- Ideal for NPOC TOC analyzers
- Performs 24/7 if needed Wall-mountable for small labs
- Few moving parts means little maintenance required and ensures long life of the generator
- Completely silent in operation
- Minimal set-up required
- 12 month comprehensive on-site warranty

## Technical Specifications

### TOC 1500 Total Organic Carbon Generator

Gas Type	Air
Max Gas Flow:	1500cc/min
Max Output Pressure:	80psi/5.5bar
Max Output Pressure:	80psi/5.5bar
Max Purity:	<1 ppm CO <sub>2</sub> / 2.0ppm THC
Min Inlet Pressure:	90psi/6.2bar
Max Inlet Pressure:	145psi/10bar
Water Purity Requirements:	N/A
Water Consumption:	TBC
Voltage:	100-230 ± 10% VAC
Frequency:	50 / 60 Hz
Current:	8 Amps
Accreditations:	CE
Size (HxWxD) mm	660 x 250 x 160 mm
Size (HxWxD) Inches	26 x 9.8 x 6.3 inches
Generator Weight	12kg / 26.4lbs

# [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

