

Infinity 5060 Nitrogen Generator

Application: N₂ supply for range of laboratory applications including LC-MS


Part number: 13-5060



Description

The Infinity 50 Series compressorless nitrogen generator utilizes membrane technology to deliver higher nitrogen flows for multiple LC-MS or LC-MS/MS instruments, Sample Evaporator and NMR applications, utilising an existing in-house supply of compressed air.

- **Gas types** Nitrogen
- **Integrated compressor:** No
- **Application Types:** N₂ supply for range of laboratory applications including LC-MS



Key Features

- Delivers laboratory-grade nitrogen at up to 360 L/min
- 24/7 operation at optimum performance if required
- Gas is supplied on demand so generator works to your schedule
- Few moving parts means little maintenance required and ensures long life of the generator
- Minimal set-up required
- Completely silent in operation
- 12 month comprehensive on-site warranty

Technical Specifications

Infinity 5060 Nitrogen Generator	
Gas Type	Nitrogen
Max Gas Flow:	360L/min
Max Output Pressure:	100psi/6.9bar
Max Output Pressure:	100psi/6.9bar
Min Inlet Pressure:	120psi/8.3bar
Max Inlet Pressure:	145psi/10bar
Min Air Inlet Flow:	180
Start Up Time:	60 mins
Voltage:	230 ± 10% VAC
Frequency:	50 Hz
Current:	0.123 Amps
Max Operating Temp:	30°C / 86°F
Size (HxWxD) mm	1150 x 500 x 800 mm
Size (HxWxD) Inches	45.3 x 19.7 x 31.5 inches
Generator Weight	95kg / 209lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001



0003