

ZA180A

Part number: 30-1660 (230v) - UK | 30-1670 (230v) - EU | 30-1680 (110v) - US



Description

Generating zero air for detector gas in GC applications, the ZA180A generator removes all hydrocarbons from a compressed air source before it enters your instrumentation, ensuring the the tidiest results from your detector.




Generator Details

- **Gas types:** Air
- **Integrated compressor:** Yes

Key Features

- Produces zero air to <0.1ppm of content
- Performs 24/7 if needed
- One generator can supply multiple instruments
- Few moving parts means little maintenance required and ensures long life of the generator
- Minimal set-up required
- 12 month comprehensive on-site warranty



Technical Specifications

	ZA180A
Gas Type	Air
Max Gas Flow:	18L/min
Hydrocarbon Content:	0.1ppm
Max Purity:	0.1ppm THC
Start Up Time:	40 Minutes
Power Consumption:	770 Watts
Voltage:	110 / 230 ± 10% VAC
Frequency:	60 Hz
Current:	2/4 Amps
Max Operating Temp:	30°C / 86°F
Accreditations:	CE
Size (HxWxD) mm	880 x 430 x 410 mm
Size (HxWxD) Inches	34.6 x 16.9 x 16.1 inches
Generator Weight	64kg / 140.8lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

