

TOC 1500 Total Organic Carbon Generator

Application: TOC Analyzers


Part number: 42-1030



Description

Designed to provide a reliable air supply for TOC analyzers, the TOC1500HP is a robust system used in labs to produce clean, dry CO₂ and CH₄ free

- **Gas types** Air
- **Integrated compressor:** No
- **Application Types:** TOC Analyzers



Key Features

- Produces clean, dry, CO₂ and CH₄-free air
- Ideal for NPOC TOC analyzers
- Performs 24/7 if needed Wall-mountable for small labs
- Few moving parts means little maintenance required and ensures long life of the generator
- Completely silent in operation
- Minimal set-up required
- 12 month comprehensive on-site warranty

Technical Specifications

	TOC 1500 Total Organic Carbon Generator
Gas Type	Air
Max Gas Flow:	1500cc/min
Max Output Pressure:	80psi/5.5bar
Max Output Pressure:	80psi/5.5bar
Max Purity:	<1 ppm CO ₂ / 2.0ppm THC
Min Inlet Pressure:	90psi/6.2bar
Max Inlet Pressure:	145psi/10bar
Water Purity Requirements:	N/A
Water Consumption:	TBC
Voltage:	100-230 ± 10% VAC
Frequency:	50 / 60 Hz
Current:	8 Amps
Accreditations:	CE
Size (HxWxD) mm	660 x 250 x 160 mm
Size (HxWxD) Inches	26 x 9.8 x 6.3 inches
Generator Weight	12kg / 26.4lbs

[**PEAK** Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK

Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America

Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China

Tel: +86 21 5079 1190
Fax: +86 21 5079 1191

For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001



0003