

## Infinity 1031 nitrogen generator

Application: 1 x SCIEX LC-MS/MS


Part number: 13-1270



### Description

The Infinity 1031 delivers high purity nitrogen and dry air to support Sciex LC-MS systems. A compressor-free solution, the Infinity 1031 requires an air source to operate. With few moving parts this generator has minimal maintenance requirements and can operate 24 hours a day for hassle free nitrogen.

- **Gas types** Nitrogen, Air
- **Integrated compressor:** No
- **Application Types:** 1 x SCIEX LC-MS/MS



### Key Features

- Capable of supplying most Sciex LC-MS/MS systems
- 24/7 operation at optimum performance
- Gas is supplied on demand so generator works to your schedule
- Few moving parts means little maintenance required and ensures long life of the generator
- Minimal set-up required
- Completely silent in operation
- 12 month comprehensive on-site warranty

## Technical Specifications

	Infinity 1031 nitrogen generator
Gas Type	Nitrogen
Max Gas Flow:	19
Max Output Pressure:	110psi/7.6bar
Gas Type	Air
Max Gas Flow:	26L/min
Max Output Pressure:	110psi/7.6bar
Gas Type	Air
Max Gas Flow:	25L/min
Max Output Pressure:	110psi/7.6bar
Max Output Pressure:	110psi/7.6bar
Min Inlet Pressure:	120psi/8.3bar
Max Inlet Pressure:	145psi/10bar
Min Air Inlet Flow:	115
Min Inlet Air Quality:	1.4.1
Start Up Time:	30 mins
Voltage:	N/A
Frequency:	N/A
Current:	N/A
Max Operating Temp:	30°C / 86°F
Particles:	<0.01 µm
Accreditations:	CE
Size (HxWxD) mm	760 x 365 x 200 mm
Size (HxWxD) Inches	29.9 x 14.4 x 7.9 inches
Generator Weight	14kg / 30.8lbs

## [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001



0003