

# Infinity 9030 Nitrogen Generator


Part number: 13-9030



## Description

The Infinity 9030 compressor-less nitrogen generator utilizes membrane technology to deliver nitrogen for LC-MS, sample evaporators, turbovap and other sample concentrators.

- **Gas types** Nitrogen
- **Integrated compressor:** No



## Key Features

- Delivers laboratory-grade nitrogen at up to 420 L/min
- 24/7 operation at optimum performance if required
- Gas is supplied on demand so generator works to your schedule
- Few moving parts means little maintenance required and ensures long life of the generator
- Minimal set-up required
- Completely silent in operation
- 12 month comprehensive on-site warranty

## Technical Specifications

Infinity 9030 Nitrogen Generator	
Gas Type	Nitrogen
Max Gas Flow:	420L/min
Max Output Pressure:	100psi/6.9bar
Max Output Pressure:	100psi/6.9bar
Min Inlet Pressure:	125psi/8.6bar
Max Inlet Pressure:	145psi/10bar
Start Up Time:	60 mins
Voltage:	230 ± 10% VAC
Current:	0
Noise Level:	10
Size (HxWxD) mm	1150 x 500 x 800 mm
Size (HxWxD) Inches	45.3 x 19.7 x 31.5 inches
Generator Weight	80kg / 176lbs

## [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

### Peak Scientific UK

Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

### Peak Scientific North America

Tel: +1-800-767-6532  
Fax: +1-978-608-9503

### Peak Scientific China

Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191

For a full list of our worldwide office locations, please visit:

Web: [www.peakscientific.com](http://www.peakscientific.com)

Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001



0003