

Corona Air Compressor

Application: Air compressor unit compatible with Corona Nitrogen units


Part number: 85-0555 (230v) | 85-1555 (110v)



Description

The Corona Air Compressor has been engineered to stack neatly with the Corona 1010 Nitrogen generator providing compressed air for laboratories without an external source of gas. The Corona Series has been developed to meet the requirements of the Thermo Scientific Corona Veo Charged Aerosol Detector.

- **Gas types** Air
- **Integrated compressor:** Yes
- **Application Types:**
Air compressor unit compatible with Corona Nitrogen units



Key Features

- Stackable with Corona 1010 Nitrogen Provides air for Corona system for labs without an external supply
- Gas is supplied on demand so generator works to your schedule
12 month comprehensive on-site warranty

Technical Specifications

	Corona Air Compressor
Gas Type	Air
Max Gas Flow:	30L/min
Max Output Pressure:	110psi/7.6bar
Max Output Pressure:	110psi/7.6bar
Min Inlet Pressure:	0
Max Inlet Pressure:	110psi/7.6bar
Power Consumption:	572 VA / 736 VA
Voltage:	110 / 230 V
Frequency:	50 / 60 Hz
Current:	5.2 A / 3.2 A
Heat output:	1950 VA / 2510 VA
Max Operating Temp:	35°C / 95°F
Particles:	<0.1µm
Accreditations:	CE, FCC
Size (HxWxD) mm	406 x 417 x 540 mm
Size (HxWxD) Inches	16 x 16.4 x 21.3 inches
Generator Weight	42kg / 92.4lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001



0003