

ZA070A

Application: Specialist Solutions


Part number: UK: 30-1610 (230v) | EU: 30-1620 (230v) | US: 30-1630 (110v)



## Description

Generating zero air for detector gas in GC applications, the ZA070A generator removes all hydrocarbons from a compressed air source before it enters your instrumentation, ensuring the the tidiest results from your detector.

- **Gas types** Air
- **Integrated compressor:** No



## Key Features

- Produces a constant flow of Zero Grade Air with a hydrocarbon content of less than 0.1ppm
- Performs 24/7 if needed
- One generator can supply multiple instruments
- Few moving parts means little maintenance required and ensures long life of the generator
- Minimal set-up required
- 12 month comprehensive on-site warranty

## Technical Specifications

	ZA070A
Gas Type	Air
Max Gas Flow:	7L/min
Hydrocarbon Content:	0.1ppm
Max Purity:	0.1ppm THC
Start Up Time:	40 Minutes
Voltage:	230 ± 10% VAC
Frequency:	50 HZ
Current:	4.9 Amps
Max Operating Temp:	30°C / 86°F
Particles:	0.01µm
Size (HxWxD) mm	880 x 430 x 410 mm
Size (HxWxD) Inches	34.6 x 16.9 x 16.1 inches
Generator Weight	64kg / 140.8lbs

# [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001



0003